



AMETEK Brookfield Launches DV Ring Light™: Illuminate Precision in Every Measurement

March 25, 2026

Middleborough, MA – AMETEK Brookfield, a global leader in precision measurement instruments, proudly announces the launch of the **DV Ring Light**, a purpose-built accessory designed to enhance visibility, confidence, and consistency in viscosity testing. Engineered exclusively for Brookfield viscometers, the DV Ring Light delivers bright, uniform illumination directly at the sample, improving clarity in low-light lab environments and supporting accurate, repeatable measurements.

DV Ring Light: Purpose-Built Illumination for Brookfield Viscometers

Designed to integrate seamlessly with Brookfield's full range of current and legacy DV and analog viscometers, the DV Ring Light transforms everyday testing conditions by improving spindle visibility and reducing shadows during immersion and measurement.

Key features include:

- **Universal Compatibility***: Works with all current and legacy DV and analog Brookfield viscometers
- **Integrated USB Power with Auto On**: Simple plug-and-play operation with automatic startup (U.S. wall adapter available)
- **Adjustable Illumination**: Built-in dimmer delivers bright, uniform lighting optimized for real lab conditions
- **Slim, Ergonomic Design**: Compact profile with intuitive left-side controls for easy access during testing
- **Instrument Protection**: Screw protectors help prevent damage to the viscometer pivot cup

*Cannot be mounted in conjunction with other pivot-cup-mounted accessories.

Why Choose the DV Ring Light?

The DV Ring Light is designed for laboratories that value precision, efficiency, and user-centric design. Unlike generic lab lighting solutions, it is engineered specifically for Brookfield instruments, ensuring optimal illumination without interfering with workflows, methods, or existing equipment. By improving visibility at the point of measurement, the DV Ring Light helps users work more confidently, supports photo-quality documentation, and enhances overall testing consistency – especially in challenging or low-light environments.

“The DV Ring Light is a simple but impactful enhancement to the viscosity testing experience,” said Chris Murray, AMETEK Brookfield. “It’s an example of how we continue to listen to our customers and design accessories that make real, practical improvements in day-to-day lab work.”

Availability

The **DV Ring Light** is available to order now and is fully supported as a genuine AMETEK Brookfield accessory. Customers are encouraged to include the DV Ring Light with new viscometer purchases or upgrade existing instruments to improve visibility and usability with minimal investment.

For more information, visit www.brookfieldengineering.com or call 508-946-6200.

About AMETEK Brookfield

With a legacy dating back to 1934, AMETEK Brookfield is synonymous with viscosity as the pioneer of the viscosity measurement instrument that set the global standard still recognized today. Since then, we have been at the forefront of innovative solutions to offer an extensive range of scientific lab instrumentation, including viscometers, rheometers, powder flow testers, texture analyzers, Computrac® moisture analysis instruments, and Jerome® toxic gas analyzers that consistently meet regulatory regulations and customer application needs. Our commitment to quality, accuracy, and customer support has made us a recognized global leader. When you choose AMETEK Brookfield, you can be confident in our instrument's repeatability, reliability, and accuracy, empowering you to make informed decisions and achieve optimal results.



Repeatable



Reliable



Accurate